

■ -80 °C Automatic Biobanking System



■ Key Features

- ✓ Refrigeration: The fully automatic sample library adopts the new double-machine overlapping refrigeration technology, which has the characteristics of low power consumption and low noise.
- ✓ Library: Intelligent library design, with the characteristics of high floor area ratio and modularization, can flexibly match various sites.
- ✓ Capacity: The storage capacity of 200,000 to 10 million 0.75ml vials.
- ✓ Versatility: Compatible storage of mainstream frozen storage vials and cryobox.
- ✓ Safety: Two sets of independent refrigeration units (one standby and one used), and equipped with liquid nitrogen refrigeration system.
- ✓ Sample management system: highly intelligent, information, support HIS\LIMS, etc.

■ Product Parameter

Product model	Biobanking sizes	Vials sizes	Sample capacity (0.75ml)
RBS50	4.6*3.2*2.9 (Customized sizes available)	0.75mL	50 million
RBS100	4.6*4.3*2.8 (Customized sizes available)	0.75mL	100 million
RBS200	11.8*3.2*2.9 (Customized sizes available)	0.75mL	200 million



Tofflon Tofflon Science and Technology Group Co.,Ltd.

Address: No.1509,Duhui Road,Shanghai,China 201108

Tel: +86 21 6490 1123 / 6490 6201

Fax: +86 21 6490 5148 / 6490 6202

E-mail: lifescience.info@tofflon.com

www.tofflon.com